DECLARACIÓN DE CONFORMIDAD CE  
CE DECLARATION OF CONFORMITY  
DÉCLARATION DE CONFORMITÉ CE

El fabricante | The manufacturer | Le fournisseur:

AUXILAB, S.L.

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*This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.*

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## 6. TROUBLESHOOTING

TROUBLE	CAUSE	SOLUTION
After turning the equipment on the display does not respond.	Power failure	Check the proper connection to power supply
	Fuse blown out	Change the fuse
	Switch failure	Contact Nahita Technical Assistance Department
Evident difference between real and display temperature	Temperature out of calibration	Calibrate temperature
	Sensor or TE module connection failure	Contact Nahita Technical Assistance Department
"OPEN" is displayed and alarm sounds	Sensor is disconnected	Contact Nahita Technical Assistance Department
"SHORT" is displayed and alarm sounds	Sensor failure	Contact Nahita Technical Assistance Department
No heating or cooling of the block	Switch power failure	Contact Nahita Technical Assistance Department
Buttons do not respond	Keyboard failure	Contact Nahita Technical Assistance Department
Block temperature decreases very slowly or it can not fall below room temperature	Fan failure	Contact Nahita Technical Assistance Department

earth connection and the socket should be handy and ready to unplug the equipment in case of emergency.

◆ Never try to repair the equipment by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the thermoblock (burns, hurts...) or damages in nearby equipments.

◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

### Cleaning

◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the thermoblock and produce an early ageing of the equipment.

◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

**ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.**

### INSTRUCTIONS ON ENVIRONMENTAL PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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


**Note:** to cancel the calibration process, press simultaneously both buttons (▲▼) for segment selection (2.4).

If you are not using the equipment for a long period of time please make sure it is disconnected from the net and protected from dust (this way you will avoid accidents and will extend its working-life).

### Seguridad

- ◆ The thermoblock must be used by previously qualified staff that knows how the equipment works thanks to the user manual.
- ◆ You should put the thermoblock in a horizontal plane stable table, having a safety area of at least 30 cm per side. Do not place the equipment near any warm supply (burners, blowlamps, etc), nor expose it directly to the sun. Avoid vibrations, dust and dry environments.
- ◆ During operation dangerous materials such as flammable or pathological substances must be out of the working area.
- ◆ When you are not using the thermoblock for a long period of time please make sure it is unplugged in order to avoid possible accidents.
- ◆ It is essential to have the equipment switched off and unplugged from the net before cleaning, checking components or replacing any piece (e.g. replacement of a fuse).
- ◆ Never try to repair the thermoblock by yourself, since you will lose the warranty and may provoke damages to the general operating system or the electrical installation, as well as injuries to the people that usually handle the equipment (burns, hurts...).
- ◆ Try not to spill any liquid on the control panel, though it is properly insulated. In case you have any doubts do please immediately unplug the thermoblock (see Maintenance).
- ◆ Made under the European regulations for electrical security, electromagnetic compatibility and security on machines.

## 5. MAINTENANCE AND CLEANING

To get the best results and a higher duration of this equipment it is essential to follow the processes of use.

**Note:** All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ The thermoblock is supplied with a Schuko standard wire. It has to be plugged to an

◆ Wait until the temperature of the thermometer remains constant for at least 20 min. After this time, if for example the thermometer marks 38.8 °C, use the temperature adjustment buttons (2.5) to adjust the displayed temperature to that shown by the thermometer. Press Start/Stop (2.3) to confirm temperature correction.

P: 40.1	ADJ *
AdjTemp=	38.8

*Note: temperature correction must be performed when temperature has remained constant for at least 20 min.*

◆ The equipment automatically enters the second calibration point; the real block temperature (P) and the second calibration point temperature (100 °C) will be displayed.

P: 40.1	ADJ *
AdjTemp=	100.0

◆ The real block temperature (P) will automatically increase until reaching 100 °C and the symbol “\*” will start flickering at the upper right side of display. When the block reaches the second calibration point temperature the symbols “ADJ” and “\*” will flicker at the same time.

P: 100.1	ADJ *
AdjTemp=	100.0

◆ Wait until the temperature of the thermometer remains constant for at least 20 min. After this time, if for example the thermometer marks 98.6 °C, use the temperature adjustment buttons (2.5) to adjust the displayed temperature to that shown by the thermometer.

P: 100.1	ADJ *
AdjTemp=	98.6

*Note: temperature correction must be performed when temperature has remained constant for at least 20 min.*

◆ Press Start/Stop (2.3) to confirm temperature correction and go out from the calibration mode.

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*Note: the protocol will always begin with segment S1 (this can not be changed) and will carry out consecutively the different segments until the last one selected. For example: if S3 is the last segment selected, the equipment will carry out segments S1, S2 and S3 consecutively.*

- ◆ Press Start/Stop (2.3) to begin the selected programmed; once finished an acoustic signal will sound.

*Note: to stop operation before the selected program has passed, press Start/Stop (2.3) for a few seconds.*

### Temperature calibration

The temperature of the equipment is already calibrated, but if desviations between real and displayed temperature are observed it can be calibrated by user.

*Note: to assure a good precision, the equipment uses a lineal temperature calibration method in two points: 40 °C and 100 °C. The temperature precision will be  $\pm$  0.5 °C after calibration.*

*Important!: during calibration both room and block temperature must be lower than 35 °C.*

Proceed as follows to perform calibration:

- ◆ Turn the equipment on and make sure that real temperature of the block is lower than 35 °C. Otherwise, let it cool.
- ◆ Fill one of the central wells of the block with olefin oil and then put a thermometer inside the well (make sure that the thermometer precision is 0.1 °C and that the bulb is completely submerged in oil). The block must be conveniently isolated so as it is not affected by environmental conditions.
- ◆ Press simultaneously both buttons (▲▼) for segment selection (2.4). The real block temperature (P) and the first calibration point (40 °C) will be displayed.

P: 20.5	ADJ
AdjTemp=	40.0

- ◆ The real block temperature (P) will automatically increase until reaching 40 °C and the symbol “\*” will start flickering at the upper right side of display. When the block reaches the first calibration point temperature the symbols “ADJ” and “\*” will flicker at the same time.

P: 40.1	ADJ *
AdjTemp=	40.0

- ◆ By pressing the corresponding buttons (2.4) select the working segment (from S1 to S5) to be programmed.
- ◆ Then adjust the working temperature value within a range from -10 °C to 100 °C by pressing the adjustment buttons (2.5) on the control panel.

*Note: when pressing temperature adjustment buttons (2.5), the cursor will automatically be placed under temperature parameter allowing its adjustment.*

- ◆ Then adjust the working time value within a range from 1 min to 99 h 59 min by pressing the adjustment buttons (2.6) on the control panel.

*Note: when pressing time adjustment buttons (2.6), the cursor will automatically be placed under time parameter allowing its adjustment.*

- ◆ Press Start/Stop (2.3) to begin operation; when real temperature of block reaches set temperature count down of programmed time will start (the symbol “:” will start flickering).

At any moment, the working segment and the programmed time and temperature (parameters at the low) as well as real temperature and remaining time (parameters at the top) will be displayed.

S3	36.8	01:24
S3	37.0	02:00

Once programmed time has passed, an acoustic signal will sound.

*Note: to stop operation before the programmed time has passed, press Start/Stop (2.3) for few seconds.*

To carry out protocols composed of several consecutive working segments, proceed as follows:

- ◆ Select time and temperature of the different working segments according to previously indicated.
- ◆ Press Prog (2.2); the following will be displayed:

Program to validate
Start: S1 End: S4

- ◆ By pressing the segment selection buttons (2.4), select the last segment (End) of the protocol to be carried out.

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- Net wire
- Allen key
- User's manual
- Warranty certificate

**We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.**

#### Installation

Before using this instrument, it is convenient for you to familiarize with its components and basic essentials.

**TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.**

♦ Please put the thermoblock on top of a horizontal, plane and stable table making a free space at least at 30 cm per side. Do not put the thermoblock near any warm supply (burners, blowlamps...), nor expose it directly to the sun, etc.

♦ Avoid the presence of flammable or toxic substances in the working area.

♦ Put the metallic block; insert the 4 screws and tighten them with the supplied Allen key.

♦ The thermoblock is supplied with a Schuko standard wire. Please insert the wire that feeds the AC electric current in the base of current 220 V, 50 Hz  $\pm 10\%$  provided with earth wire and to the other end to the thermoblock connector (1.4).

**Neither the manufacturer nor the distributor will assume any responsibility for the damages produced to the equipment during its installation or damages to people suffered by the improper use of the electric connection. The tension should be 220 V, 50 Hz  $\pm 10\%$ .**

#### Operation mode

The equipment can operate in working segments (up to 5 segments) each one of which allows setting time (from 1 min to 99 h 59 min) and temperature (from -10 °C to 100 °C); moreover, the equipment can be programmed to carry out protocols composed of several consecutive segments.

♦ Turn the equipment on by pressing power on/off (1.3). After a few seconds the following parameters (or similar) will be displayed; the left column makes reference to the working segment, the middle to the temperature (°C) and the right to the time (h:m).

S1	15.0	00:00
S1	37.0	02:00

### 3. TECHNICAL SPECIFICATIONS

Code	50603020
Model	603/20
Temperature range *	-10 °C – 100 °C
Temperature accuracy	≤ 0.3 °C
Temperature resolution	0.1 °C
Heating time (from 20 °C to 100 °C)	≤ 20 min
Cooling time (from 20 °C to -10 °C)	≤ 30 min (if RT is 25 °C)
Timer	1 min – 99 h 59 min
Working condition	5 °C – 30 °C; ≤ 70% HR
Consumption	150 W
Power	CA 220 V, 50-60 Hz
Dimensions (LxWxH)	30x21x18 cm
Weight	5 Kg

\* The thermoblock can reach a temperature of -10 °C if room temperature is lower than 25 °C.

The thermoblock is supplied without metallic block so as each user can make up the equipment according to their specific necessities. The following blocks are available:

Code	Capacity
90603021	96 x 0.2 mL
90603022	54 x 0.5 mL
90603023	35 x 1.5 mL
90603024	35 x 2.0 mL

### 4. INSTALLTION / SETTING UP

#### Preliminary inspection

Unwrap the thermoblock, take off the involving plastic and take off the polyspan protection in which it comes fitted. Take off all the protective items and, without connecting the equipment to the net, make sure that it does not present any damage because of the shipment. In case the thermoblock presents any damage tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

Please check that all the accessories are enclosed with the equipment:

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

## 2. DESCRIPTION


- 1.1 Metallic block (not included)
- 1.2 Control panel
- 1.3 Power on/off

- 1.4 Socket
- 1.5 Fuse

- 2.1 LCD display
- 2.2 Prog: operating program selection
- 2.3 Start/Stop
- 2.4 Seg: segment selection buttons

- 2.5 Temp: temperature adjustment buttons
- 2.6 Time: time adjustment buttons

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## 1. USES OF THE INSTRUMENT

The Nahita thermoblock model 603/20 allows incubating samples under conditions of total and precise temperature control thus, being useful in Molecular Biology and Genetics for DNA sample analysis, extraction of lipids and other cellular components, denaturalization of proteins previous to electrophoresis, PCR amplification, etc.

### Main features

- Time and temperature digital control and display in LCD.
- Real and set temperatures are displayed simultaneously.
- Possibility of using different metallic blocks to adapt to specific necessities in each protocol. Blocks are easily interchangeable and very simple to clean and sterilize.
- Temperature can be calibrated to adjust to user necessities.
- Acoustic signal at the end of the selected operation program

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita thermoblock 603/20. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

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Please bear in mind the following:

- ◆ This manual is inseparable from the Nahita thermoblock 603/20, so it should be available for all the users of this equipment.
- ◆ You should carefully handle the thermoblock avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ◆ Never dismantle the different pieces of the thermoblock to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the thermoblock.
- ◆ To prevent fire or electric discharges avoid dry or dusty environments. In case it may happen unplug the equipment immediately.
- ◆ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by email to [info@auxilab.com](mailto:info@auxilab.com) or contacting Nahita Technical Assistance Department by telephone: +34 807 117 040 (1.30 Euros/min).
- ◆ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ◆ Overhaul is not covered by the thermoblock warranty.
- ◆ Operations made by non-qualified staff will automatically produce a loss of the equipment warranty.
- ◆ Neither fuses nor accessories (including their loss), are covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ◆ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ◆ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ◆ Manufacturer reserves the right to modify or improve the manual or equipment.